



Leica CM1850

The cryostat for routine histology
and clinical pathology

Ergonomic design with outstanding safety features

For 125 years Leica has been providing state-of-the-art microtome technology and this commitment is still the driving force behind our new developments. The design of the Leica CM1850 cryostat incorporates the proven excellent features of its predecessor, the Leica CM1800, but additional emphasis was placed on power saving, increased efficiency and operator safety. The result is a versatile cryostat with optimized cooling system, rapid specimen freezing and smooth specimen orientation for the high-quality sectioning demanded in routine histology and clinical pathology. Up-to-date technology with an ergonomic design – for quick and safe sectioning.



Spacious cryochamber

The perfectly illuminated cryochamber is easily accessible through a large heated sliding window. Ample space is provided for convenient specimen handling. Smooth surfaces inside ensure safe and easy cleaning.

Enhanced insulation

A new insulation technology leads to power savings of approximately 10% compared to conventional heat insulation systems. This new feature enhances the durability of the refrigerating system and ensures stable cryochamber temperatures, even when operating in unfavorable environmental conditions in the lab such as high humidity.

Easy-to-service

In developing the Leica CM1850 emphasis was placed on reliable performance and a service-friendly design. All serviceable components are easily accessible.



Quick specimen freezing

An actively cooled freezing shelf allows extremely fast specimen freezing down to -60 °C with the Peltier unit. Manual defrosting prevents icing of the freezing shelf ensuring that it is always ready to use.



Functional CE knife holder

The lateral displacement feature of the CE knife holder for disposable blades allows the use of the entire blade length without need to reposition the blade. A glass anti-roll guide integrated into the holder is automatically positioned with blade movement. With the new CE knife holder serial sections down to 3 µm are easily obtained.



Easy handling and disinfection

The microtome is protected by a slot cover to allow easy and safe spray disinfection with Leica Cryofect. The smooth specimen orientation with automatic centering and zero notch features outstanding handling capabilities for fast and precise sectioning.

Self-explanatory control panel

The control panel is divided into functional groups. Large, easy-to-read LEDs ensure ease of operation and prevent operating errors.



Enclosed drainage system

A closed drainage system ensures that condensation and defrost fluids are collected into a container placed at the front of the cryostat cabinet allowing controlled disposal of fluids. This considerably reduces the release of germs into the air.



Ergonomic coarse feed panel

The motorized coarse feed control is ergonomically positioned in the left arm rest. Either fast or slow specimen advance and return can be selected.



At a glance...

- Chamber temperature down to -35 °C
- Actively cooled quick freeze shelf with Peltier unit for freezing down to -60 °C.
- Power-saving refrigeration system
- Splash-protected microtome
- Motorized coarse feed at two speeds
- Smooth specimen orientation with automatic centering
- Automatic and manual chamber defrost cycles
- Manual quick freeze shelf defrost cycle
- Ergonomic design with convenient foot rest
- LED
- Compact dimensions
- Efficient service and maintenance concept



Leica CM1850 – Technical specifications

Section thickness setting:	1 to 60 µm	Cooling capacity:	
	0 - 10 µm in 1 µm increments	Temperature setting range:	0 °C to -35 °C
	10 - 20 µm in 2 µm increments	Cooling down time to -35 °C:	3,5 hours
	20 - 60 µm in 5 µm increments	Chamber defrosting:	programmable
Maximum specimen size:	55 mm Ø		1 automatic defrost cycle/24 hours
Horizontal specimen feed:	25 mm		duration: max. 12 min.
Vertical specimen stroke:	59 mm		manual defrost cycle
Specimen retraction:	standard on all versions except US-version	Freeze shelf temperature:	down to -40 °C (± 2K)
Specimen orientation:	Standard delivery	Peltier unit:	down to -60 °C (+ 5K) at maximum freeze shelf temperature
Trimming:	via motorized coarse feed	Freeze shelf defrosting:	manual defrost cycle
Motorized coarse feed:	2 speeds	Cryocabinet:	
	rapid: 0.7 mm/s	Dimensions:	670 x 1165 x 695 mm
	slow: 0.2 mm/s	Weight (incl. microtome):	131 kg

As confirmed by the successful c-CSA-us certification, the Leica CM1850 has been designed and manufactured in compliance with UL, CSA and IEC requirements.

State-of-the-art development, manufacturing and quality control


procedures – certified under DIN EN ISO 9001 – ensure highest quality and reliability.

A wide range of accessories available on request.

Technical specification subject to change without prior notice.

For further information please contact your local Leica sales company

Australia:	Gladesville/NSW	Phone +61 2 9897 9700	Fax +61 2 9817 8358
Austria:	Vienna	Phone +43 1 486 80 50	Fax +43 1 486 80 50 30
Canada:	Richmond Hill/Ontario	Phone +1 905 762 2000	Fax +1 905 762 8937
Denmark:	Herlev	Phone +45 4454 0101	Fax +45 4454 0111
France:	Rueil-Malmaison Cedex	Phone +33 1 4732 8585	Fax +33 1 4732 85 86
Germany:	Bensheim	Phone +49 6251 136 0	Fax +49 6251 136 155
Italy:	Milan	Phone +39 0257 486.1	Fax +39 0257 40 3273
Japan:	Tokyo	Phone +81 3 5421 2804	Fax +81 3 5421 2894
Korea:	Seoul	Phone +82 2 514 6543	Fax +82 2 514 6548
Netherlands:	Rijswijk	Phone +31 70 4132130	Fax +31 70 4132139
People's Republic of China:	Hong Kong	Phone +852 2 564 6699	Fax +852 2 564 4163
Portugal:	Lisbon	Phone +351 1 388 9112	Fax +351 1 385 4668
Singapore:	Singapore	Phone +65 6779 7823	Fax +65 6773 0628
Spain:	Barcelona	Phone +34 93 494 9530	Fax +34 93 494 9532
Sweden:	Sollentuna	Phone +46 8 6254 545	Fax +46 8 6254 510
Switzerland:	Glattbrugg	Phone +41 1 809 34 34	Fax +41 1 809 34 44
United Kingdom:	Milton Keynes	Phone +44 1 908 246246	Fax +44 1 908 609992
USA:	Bannockburn/Illinois	Phone +1 847 405 0123	Fax +1 847 405 0164



...and more than 100 national distributors.